IMPETUX

The Force of Light

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
IMPETUX

The Force of Light

The optical manipulation platform developed by IMPETUX consists of an optical trapping unit (**Cygnium G-422**) combined with a force sensor unit (**Lunam T-40i**).

Direct measurement of **forces inside cells and tissues**.

Calibration-free force measurements.

Active micro-**rheological information**.

Compact “on-microscope” opto-mechanical design.

Compatible with different microscopies (epi-FL, DIC, TIRFF,...)

User-friendly software.

Multiple trapping capabilities for complex manipulation experiments (>10 traps).

The **optical trapping** system **Cygnium G-422** has been specifically designed to work with living cells and tissues. It is based on a patented technology that overcomes the main limitations of current systems, bringing key advantages for the measurement of forces in complex media.

The optical tweezers Cygnium G-422 is conceived to work as a satellite component of an inverted microscope. The system offers the possibility to **work with multiple traps** (>100) simultaneously synchronized (**time sharing**) by means of acousto-optic deflection technology, **automatically measuring the force** exerted by each trap **and allowing the acquisition of fluorescence images at the same time**.

Cygnium G-422 allows the possibility to determine the rheological properties of the working medium. By mean of active/passive techniques, IMPETUX technology offers the possibility of **measuring micro-rheological magnitudes inside complex mediums like cells cytoplasm** with great viscosity and abundant of elastic properties.

The main performances and technical specifications of the optical trapping unit are:

- **Ultra-fast beam steering through acousto-optic deflectors for the generation of multiple traps (>100 independent traps).**
- 2D fast trap steering (100 kHz). **Force measurements at 100 kHz can be simultaneously performed for up to 10 independent traps.**
- Accurate positioning of the individual traps.
- Ultra-stable low noise Infrared laser for biological compatibility.
- **Optimized laser power homogeneity (<1 %)** of all traps generated in the working field (typically 80 μm x 80 μm for a 60X objective magnification).
- **Laser power at sample as high as >0.5 W.**
- **Maximum typical trapping forces of 100-500 pN.**
- **Simultaneous Fluorescence imaging.**
- **Compatible with different microscopies (DIC, TIRFF...)**
- User-friendly control and acquisition software included.
- Robust high-quality aluminum enclosure of the optical unit specially designed for the maximum mechanical stability of the system.
- Dimensions of the Optical unit (L x H x W): 360 x 250 x 90 mm.

The laser unit (included in the system) and Cygnium G-422 electronic drivers are mounted into a high quality enclosure (either rack or table top enclosure).

The **force sensor** module [Lunam T-40i](#) is specially conceived to simultaneously measure the forces exerted by several optical traps generated by the optical trapping module. **Our patented and distinctive technology** measures the force as a derived magnitude of the direction of propagation of the laser beam, instead of the sample position, which

provides a direct route to the force. After an accurate calibration at factory, the deflection of the beam produced by the sample directly corresponds to the optical force exerted on the trapped particle.

This gives our technology a key advantage over the competence systems, usually based on back focal plane interferometry, as no restrictions on the sample are imposed and no previous calibration from the user side is required. **This is especially advantageous when working with living cells.**

Among the most relevant advantages of the Lunam T-40i force sensor compared to other competitors, we highlight:

- Ability to perform **direct force measurements** and work in **sample position detection** mode.
- Permanent calibration of the system made at factory: **No need of any previous calibration to start measuring** or when sample is changed.
- **Reliable force measurements inside thick samples (up to 0.5 mm).**
- **Exact force measurement independently of sample size and shape** (including non-spherical samples).
- **Precise determination of the laser power at sample.**
- **Possibility of measuring beyond the linear region of the trap:** under the same conditions, Lunam T-40i is capable of measuring higher forces (a 2-3 fold factor) than its competitors. This **significantly reduces the laser power irradiation of the sample** necessary to measure a certain force.

The main performances and technical specifications of the Lunam T-40i force sensor unit are:


- **Force resolution <50 fN.**
- **Position resolution <1nm.**
- **Compatible with commercial micro-incubation chambers and confocal imaging chambers** (like for instance the RC-30WA model from Warner Instruments).
- High Numerical-Aperture (NA=1.4) immersion optics.
- Integrated sensor noise over the whole bandwidth <0.1 pN (typ.).
- Temperature-compensated, duo-lateral position sensitive detector (PSD) (max. sampling frequency 100 kHz).
- Up to 100 kHz, 18-bit, analog-to-digital conversion.
- Highly-regulated, low-noise linear power supply (models with 100/120/220/240 VAC- 50/60 Hz available).
- Symmetric and robust opto-mechanical design.
- Sensor head unit dimensions: 22 cm x 18 cm x 11 cm.
- Direct PC communication through Hi-speed USB 2.0 port.
- User-friendly acquisition software and LabVIEW libraries included.
- Straightforward installation and tuning routines ensure correct measurements and reproducibility.

Additionally, the Lunam T-40i force sensor can also be complemented with additional options like:

- **Force clamping** modulus to apply constant force levels on samples
- **Z axis force measurement** modulus for **3D force measures**.

The control and acquisition software supplied by Impetux consists of a **Labview application** with moduli for the **control of several microscope cameras** (like for instance SCMOS Hamamatsu ORCA Flash 4.0 V3), the generation of optical traps and the acquisition of synchronized force measurements. The **generation of the traps** as well as their positioning and displacement over the working field can be **controlled by simply clicking and dragging the mouse** over the camera screen interface. The **traps can also be generated through predefined files** in order to execute finely controlled trap displacements. The direct force measurements of each tarp generate useful data (up to 100 kHz sampling frequency) in real time which can be stored in several format files for their ulterior analysis. **Automated routines are also implemented in the software in order to perform active micro-rheological studies of the samples.**

LabVIEW is the preferred software development tool used in scientific instrumentation environments and has the advantage that it allows implementing simple modifications by the user, Impetux technicians or third-party developers in case of need. For this purpose, **IMPETUX distributes .vi modules with open source code.**

An additional advantage of LabVIEW is the ease of implementation into the software of control routines for any new equipment that is added later to the system (such as new cameras or motorized platforms) since most of the instrumentation manufacturers facilitate the control drivers of their equipment in this environment.