



# 一次性注射器和套件

提供各种1mL至140mL的一次性塑料注射器。您可以获得盒装数量的折扣。

Description (mL=CC)	Part #	Price Each	Box Part #	Box Price
140 mL Luer Lock Non-Inert	SYR-MJ140LL	询价	SYR-MJ140LLBX	询价(每盒20个)
60 mL Luer Lock Non-Inert	不再提供	询价	不再提供	询价 (每盒25个)
50 mL Luer Lock Inert	SYR-NJ50LL	询价	SYR-NJ50LLBX	询价 (每盒30个)
30 mL Luer Lock Inert	SYR-NJ30LL	询价	SYR-NJ30LLBX	询价(每盒50个)
20 mL Luer Lock Inert	SYR-NJ20LL	询价	SYR-NJ20LLBX	询价(每盒100个)
10 mL Luer Lock Inert	SYR-NJ10LL	询价	SYR-NJ10LLBX	询价(每盒100个)
5 mL Luer Lock Inert	SYR-NJ5LL	询价	SYR-NJ5LLBX	询价 (每盒100个)
3 mL Luer Lock Inert	SYR-NJ3LL	询价	SYR-NJ3LLBX	询价 (每盒100个)
1 mL Luer-Lock Non-Inert	SYR-EX1LL	询价	SYR-EX1LLBX	询价(每盒100个)
1 mL Tuberculin Inert	SYR-TB1	询价	SYR-EX1LLBX	询价 (每盒100个)
Large Syringe kit	P-SYRKIT-LG	询价	N/A	N/A
Small Syringe kit	P-SYRKIT-SM	询价	N/A	N/A
Dual Pump Plumbing Kit	P-DKIT	询价	N/A	N/A



Large Syringe Kit  
Part ID: P-SYRKIT-LG

- (2) 20mL Luer Lock Inert
- (2) 30mL Luer Lock Inert
- (2) 60mL Terumo Luer Lock Non-Inert
- (6) Female Luer Barbs
- (3) Male Luer Barbs
- (20 ft) Tubing .125" ID; .188" OD, PVC



Small Syringe Kit  
Part ID: P-SYRKIT-SM

- (2) 1mL Luer Lock Inert
- (2) 3mL Luer Lock Inert
- (2) 5mL Luer Lock Inert
- (2) 10mL Luer Lock Inert
- (6) Female Luer Barbs
- (3) Male Luer Barbs
- (20 ft) Tubing .125" ID; .188" OD, PVC



Dual Pump Plumbing Kit  
Part ID: P-DKIT

- For use with our continuous infusion system or dual infusion system (Dual-NE-1000). Provides a dual check valve for two syringes, withdraw will pull fluid from a reservoir and infusion will output into a separate line
- The kit is pre-assembled, all you need to do is plug in the syringes
- Each dual check valve uses 6" of 1/8" tubing to connect to a Y connector, one for reservoir and one for output; each Y-connector has 5' of tubing to make connections to your application. Two 60cc Terumo syringes are included

- Highly corrosive chemicals are not compatible with this kit. Contact us for a solution that is compatible.