

## **FLOW RATES**

The Model BSP-99M can accommodate from a 50cc syringe to a microliter syringe, both in plastic and in glass. The Model BSP-99M has a built in flow chart, calibrated for the Multifit brand syringe, with 99 equally spaced speeds from

|                         |                   |
|-------------------------|-------------------|
| .145 to 143. cc/hr..... | with 50cc syringe |
| 0.95 to 94. cc/hr.....  | with 30cc syringe |
| 0.71 to 70. cc/hr.....  | with 20cc syringe |
| 0.40 to 39. cc/hr.....  | with 10cc syringe |
| 0.25 to 25. cc/hr.....  | with 5cc syringe  |

To determine flow rates with other syringes, a calibration sheet is included with every pump.