

## FLEXCELL® FLEXLINK® *Tension or Compression*

(cat. nos. FX-5000TFL or FX-5000CFL)

### Run multiple regimens at one time!

- Additional strain unit module that links with a FX-5000 system computer allowing the user to simultaneously operate two different regimens on two different baseplates.
- Connects in a series to operate up to four FlexLinks from a central computer.

\* FlexLink may require an additional vacuum or positive pressure source depending on the intensity of the regimens.



**Figure 11.** Image of FlexLink Tension System controller.



#### **FX-5000 FlexLink includes:**

- FX-5000 FlexLink
- Accessory kit